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
**READ THIS MANUAL COMPLETELY AND
CAREFULLY BEFORE USING THIS PRODUCT**

Keep this manual in a safe location for future reference

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1. IMPORTANT SAFETY NOTES

Signs and symbols

The following signs appear in the Safety Section (pages 2, 3) and in this manual on pages 5, 6, 7.

READ AND UNDERSTAND THIS ENTIRE MANUAL INCLUDING THE SAFETY SECTION AND ALL INSTRUCTIONS AND WARNINGS COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT. FOLLOW ALL SAFETY INSTRUCTIONS AND WARNINGS TO AVOID HAZARDOUS SITUATIONS AND TO MAKE CORRECT USE OF THIS PRODUCT.

WARNING

WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

CAUTION indicates a hazardous situation which, if not avoided, may result in minor or moderate injury.

NOTICE

NOTICE addresses practices not related to personal injury, such as product and/or property damage.

WARNING:

- Packaging materials are a deadly hazard for children and can cause suffocation. Remove all packaging materials immediately and keep them away from children at all times.
- Swallowing batteries and/or battery fluid can be extremely dangerous. Keep the battery and scale out of the reach of children and disabled persons. Should any person swallow a battery and/or battery fluid, please call 911 immediately.
- Should battery fluid leak and come into contact with your eyes or skin, immediately rinse with plenty of clean water. Please call 911 immediately.
- Batteries should not be charged or reactivated by any other means.
- Batteries should not be taken apart, thrown in the fire or short circuited.
- Do not step on the edge of the scale (e.g., when getting on or off), otherwise it may tip.
- Do not jump on the scale.
- Protect the scale against hard knocks, temperature fluctuations and heat sources which are too close (e.g., stoves, heating radiators).
- Do not drop the scale or drop any objects on it.
- Do not step or stand on the scale when your body and/or feet are wet, for example after taking a bath or shower. There is a danger of slipping.
- Do not step on the scale when its surface is wet. There is a danger of slipping.
- Do not use the unit on slippery surfaces such as wet floors. There is a danger of slipping.
- This product does not and is not intended to provide a medical diagnosis. For medical advice please contact a licensed physician. If you have an emergency please call 911 immediately.
- Consult your physician or healthcare provider before beginning a weight reduction or exercise program.
- Always consult your physician to determine what will be suitable for you.
- Make sure that you enter the desired unit of mass regarding the weight results.
- This product is not a toy. Keep it out of reach of children.
- People with disabilities, such as impairments, activity limitations, and participation restrictions, or who are physically frail, should be assisted by another person when using this scale.
- Do not step on the scale if you weigh more than 400 lb (180 kg/28 st). When in doubt contact your physician.
- Do not weigh more than one person on the scale at one time.
- Do not use this product for any purpose other than weighing human beings.
- Place the scale on level flooring, otherwise it may tip.
- The product may have sharp edges. You could get injured.
- If the glass platform is damaged, there could be risk of danger or injuries caused by broken glass.

CAUTION:

- This unit is intended for home use only. It is not intended for professional use in hospitals or other medical facilities. This unit does not support the standards required for professional use.
- Please note that technically related measuring tolerances are possible, as the scale is not calibrated for professional medical applications.
- If a battery has leaked do not touch the battery fluid. Avoid skin contact (e.g., put on protective gloves) and clean the battery compartment with a dry cloth.
- Carefully step on the scale surface. Pay attention to avoid stepping with your toes under the scale surface. This might lead to injuries.
- Do not disassemble this unit. You could get injured.
- Do not put your fingers in any slot, inside or under the scale.

NOTICE:

- Protect the scale against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (e.g., stoves, heating radiators). This may damage the unit. Please see section 6 "Cleaning, maintenance and disposal" and section 8 "Technical specifications" for details.
- Do not use any aggressive solvents, cleaning agents, detergents or any other strong chemicals to clean the device.
- Never immerse and or spill water or any other liquid onto the scale, otherwise liquid will enter the scale and cause damage.
- Do not wash the scale in a washing machine, dish washer or tumble dryer.
- Always use the size and type of battery indicated.
- Never use rechargeable batteries. This may damage the unit.
- Replace weak battery before it discharges completely. When "Lo" appears in the display please change the battery.
- Leaking batteries may damage the device. If you do not intend to use the scale for longer periods, remove the battery from the battery compartment before placing the device in storage.
- To avoid draining the battery do not place any items on the scale when not in use. (Please make sure that there is no weight on the feet of the scale.)
- Batteries can contain toxins that are harmful to the environment. Always dispose of batteries in accordance with applicable local regulations. Do not dispose of batteries with normal household waste.
- Observe the local regulations for material disposal. Dispose of the device in accordance with local regulations.
- Do not step on the scale if you weigh more than 400 lb (180 kg/28 st). This may damage the unit. When in doubt contact your physician.
- Do not weigh more than one person on the scale at one time. This may damage the unit.
- Do not use mobile telephones, microwaves, and other devices which generate strong electrical or electromagnetic fields near the unit. This may result in incorrect operation of the unit.
- Never attempt to repair the instrument or adjust it yourself. This may damage the unit and impair the functions. If you need to have the unit repaired, please contact our customer service. Please see warranty for service contact. Before submitting any complaint, first check the battery and replace it if necessary.
- Changes or modifications to the device will nullify the user warranty.

Content

1. Important Safety Notes	2	4. Batteries	5
▲ WARNING	2	5. Operation	5
▲ CAUTION	3	6. Cleaning, maintenance and disposal	6
NOTICE	3	Cleaning and maintaining	7
2. Package contents	4	Disposal	7
3. Interesting facts	4	7. FCC statement	7
3.1 Weighing correctly	4	8. Technical specifications	7
3.2 General notes	4	9. Troubleshooting	8
		10. Warranty	8

SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE

If the instruction manual is damaged or if you no longer have the instruction manual in your possession, please contact customer service at 1-800-536-0366.

Dear Customer,

Congratulations on your purchase of the CVS/pharmacy® Glass Scale.

Please read these instructions completely and carefully before using this product and keep them for later reference. Be sure to make them accessible to other users and observe the information they contain.

2. Package contents

- GS05 Glass Scale
- 1 x 3 V Lithium Battery CR2032
- Operating Instructions

3. Interesting facts

3.1 Weighing correctly

Place the scale on a hard, plane level floor (no carpet); a firm floor covering is required for correct measurement.

- If possible, always weigh yourself at the same time of day (preferably in the morning) striving for comparable conditions such as after using the toilet, on an empty stomach and unclothed, in order to achieve comparable results.
- Stand still and upright during measurement.

3.2 General notes

- The capacity of the scale is 400 lb (180 kg/28 st).
- The results are shown in 0.2 lb increments (100 g/0.02 st).
- On delivery the scale is set to the unit "lb." A switch located on the back of the scale enables you to change to "kg" or "st."

Measuring units

lb	pound
kg	kilogram: the base unit of mass in the International System of Units
st	stone: only for countries where British English is spoken

For product information, questions or customer service call toll free 1-800-536-0366.

4. Batteries

WARNING:

- Swallowing batteries and/or battery fluid can be extremely dangerous. Keep the battery and scale out of the reach of children and disabled persons. Should any person swallow a battery and/or battery fluid, please call 911 immediately.
- Should battery fluid leak and come into contact with your eyes or skin, immediately rinse with plenty of clean water. Please call 911 immediately.
- Batteries should not be charged or reactivated by any other means.
- Batteries should not be taken apart, thrown in the fire or short circuited.


CAUTION:

- If a battery has leaked do not touch the battery fluid. Avoid skin contact (e.g., put on protective gloves) and clean the battery compartment with a dry cloth.

NOTICE:

- Always use the size and type of battery indicated.
- Never use rechargeable batteries. This may damage the unit.
- Replace weak battery before it discharges completely. When “Lo” appears in the display please change the battery.
- To avoid draining the battery do not place any items on the scale when not in use. (Please make sure that there is no weight on the feet of the scale.)
- Leaking batteries may damage the device. If you do not intend to use the scale for longer periods, remove the battery from the battery compartment before placing the device in storage.
- Batteries can contain toxins that are harmful to health and the environment. Always dispose of batteries in accordance with applicable local regulations. Do not dispose of batteries with normal household waste.

If present, remove the battery insulating strip or remove the battery's protective film and insert the battery according to the polarity. If the scale fails to operate, remove the battery completely and reinsert it. Replace the battery cover carefully.

A new battery will last for approximately 1.5 years if the unit is used three times per day.

The trial battery supplied with the unit may have a shorter life. Remove battery if you do not intend to use the scale for longer periods before placing the device in storage.

Your scale is equipped with a “replace battery” indicator. When the scale is operated with a battery that is too weak, “Lo” appears in the display field and the scale automatically switches off. The battery must be replaced in this case (1 x 3 V lithium battery CR2032).

5. Operation

WARNING:

- Packaging materials are a deadly hazard for children and can cause suffocation. Remove all packaging materials immediately and keep them away from children at all times.
- Do not step on the edge of the scale (e.g., when getting on or off), otherwise it may tip.
- Do not jump on the scale.
- Protect the scale against hard knocks, temperature fluctuations and heat sources which are too close (e.g., stoves, heating radiators).
- Do not drop the scale or drop any objects on it.
- Do not step or stand on the scale when your body and/or feet are wet, for example after taking a bath or shower. There is a danger of slipping.
- Do not step on the scale when its surface is wet. There is a danger of slipping.
- Do not use the unit on slippery surfaces such as wet floors. There is a danger of slipping.
- This product does not and is not intended to provide a medical diagnosis. For medical advice please contact a licensed physician. If you have an emergency please call 911 immediately.
- Consult your physician or healthcare provider before beginning a weight reduction or exercise program.
- Always consult your physician to determine what will be suitable for you.

- Make sure that you enter the desired unit of mass regarding the weight results.
- People with disabilities, such as impairments, activity limitations, and participation restrictions, or who are physically frail, should be assisted by another person when using this scale.
- Do not step on the scale if you weigh more than 400 lb (180 kg/28 st). When in doubt contact your physician.
- Do not weigh more than one person on the scale at one time.
- Do not use this product for any purpose other than weighing human beings.
- Place the scale on level flooring, otherwise it may tip.
- The product may have sharp edges. You could get injured.
- If the glass platform is damaged, there could be risk of danger or injuries caused by broken glass.

CAUTION:

- This unit is intended for home use only. It is not intended for professional use in hospitals or other medical facilities. This unit does not support the standards required for professional use.
- Please note that technically related measuring tolerances are possible, as the scale is not calibrated for professional medical applications.
- Carefully step on the scale surface. Pay attention to avoid stepping with your toes under the scale surface. This might lead to injuries.
- Do not disassemble this unit. You could get injured.
- Do not put your fingers in any slot, inside or under the scale.

NOTICE:

- Do not step on the scale if you weigh more than 400 lb (180 kg/28 st). This may damage the unit. When in doubt contact your physician.
- Do not weigh more than one person on the scale at one time. This may damage the unit.
- Do not use mobile telephones, microwaves, and other devices which generate strong electrical or electromagnetic fields near the unit. This may result in incorrect operation of the unit.


This scale is used to determine the weight of human beings. The scale has a weighing range up to 400 lb (180 kg/28 st).

Place the scale on a hard, plane level floor (no carpet); a firm floor covering is required for correct measurement.

Note: For first use and relocation

In order to ensure accurate results it is recommended to reset the scale to "0."


- Step onto the scale for one second and then step off.
- The scale shows "0.0 lb" (kg/st).
- It will shut off automatically after several seconds.

Quickstart: Step onto the scale.

Stand on the scale without moving and distribute your weight evenly on both legs.

After the weight is displayed step off the scale. The scale will automatically turn off after several seconds.

6. Cleaning, maintenance and disposal

NOTICE:

- Protect the scale against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (e.g., stoves, heating radiators). This may damage the unit. Please see section 8 "Technical specifications" for details.
- Do not use any aggressive solvents, cleaning agents, detergents or any other strong chemicals to clean the device.
- Never immerse and or spill water or any other liquid onto the scale, otherwise liquid will enter the scale and cause damage.
- Do not wash the scale in a washing machine, dish washer or tumble dryer.

- Observe the local regulations for material disposal. Dispose of the device in accordance with local regulations.
- Never attempt to repair the instrument or adjust it yourself. This may damage the unit and impair the functions. If you need to have the unit repaired, please contact our customer service. Please see warranty for service contact. Before submitting any complaint, first check the battery and replace it if necessary.
- Changes or modifications to the device will nullify the user warranty.

Cleaning and maintaining

To keep your scale in the best condition and protect the unit from damage follow the directions below: The scale should occasionally be cleaned with a damp cloth. Do not use aggressive cleaning agents. The unit is not waterproof, do not immerse the unit in water and do not spill water over the unit. If your scale accidentally gets wet immediately clean off the water. Store the unit in a safe and dry location. Do not expose the unit to direct sunlight, extreme hot or cold temperatures or humidity. Thoroughly dry any moisture off the unit before storing. Store the unit safely out of the reach of children. Do not subject the unit to strong shocks, such as dropping the unit on the floor. Always use the unit following the instructions provided in this manual.

Disposal

NOTICE:

- Batteries should be disposed of separately from household waste. Always dispose of batteries as per your local regulations. Again, do not dispose of batteries with normal household waste.
- Observe the local regulations for material disposal. Dispose of the device in accordance with local regulations. If you have any questions, please contact the local authorities responsible for waste disposal.

7. FCC statement

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules.

These limits are designed to provide reasonable protection against harmful interference in a residential installation. The equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and the receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

8. Technical specifications

Battery type	1 x 3 V lithium battery CR2032
Measuring range	11 lb to 400 lb (5 kg to 180 kg / 0.79 st to 28 st)
Weight unit	0.2 lb increments (100 g/0.02 st)
Temperature operating range	41 °F to 95 °F (5 °C to 35 °C)
Storage temperature range	14 °F to 140 °F (-10 °C to 60 °C)
Scale size	11.8" x 11.8" x 0.9" (300 mm x 300 mm x 23.5 mm)
Net weight of scale	2.87 lb (1.3 kg)

Subject to technical modifications.

9. Troubleshooting

The following error messages can appear.

Error display	Possible causes of errors	Correction of errors
OL	The maximum capacity of 400 lb (180 kg/28 st) has been exceeded.	Weigh only the maximum permissible weight.
Lo	Low power.	Change the battery.
E	If the scale has not been in use for a long period of time.	Repeat measurement.
	The scale was switched on by hand but not on the floor.	Wait until the scale switches off automatically. Place the scale on level flooring (no carpet). Repeat measurement.
	No display.	Battery is empty, battery was inserted incorrectly or the scale is broken.
	Scale shows unrealistic result.	Wait until scale switches off automatically. Repeat measurement. Place the scale on level flooring (no carpet). Don't move during weighing.

10. Warranty

Limited Lifetime Warranty For Original Purchaser

Your CVS/pharmacy® Glass Scale, Model GS05, excluding the battery, is warranted to be free from defects in materials and workmanship for the life of the product under normal conditions of intended use and service. This warranty extends only to the original retail purchaser and does not extend to retailers or subsequent owners.

We will, at our option, repair or replace the CVS/pharmacy® Glass Scale, Model GS05, without additional charge, for any part or parts covered by these written warranties. No refunds will be given. Repair or replacement is our only responsibility and your only remedy under this written warranty. If replacement parts for defective materials are not available, the warrantor reserves the right to make product substitutions in lieu of repair or replacement.

For warranty service contact our customer service department at 1-800-536-0366 to provide a description of the problem. If the problem is deemed to be within the scope of the limited lifetime warranty, you will be asked to mail the product at your costs in its original package with proof of purchase, your name, address and phone number. If the problem is not deemed to be within the scope of the limited lifetime warranty, we will provide a quotation for repair respectively replacement and return shipping fee.

This warranty does not cover damage caused by misuse or abuse; accident; the attachment of unauthorized accessory; alteration to the product; improper installation; misapplication; lack of reasonable care with respect to the product; unauthorized repairs or modifications; improper use of electrical/power supply; old worn batteries; normal wear; loss of power; dropped product; malfunction or damage of an operating part as a result of failure to comply with instructions for use or to provide manufacturer's recommended maintenance; transit damage; theft; neglect; vandalism; environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of the warrantor. This warranty is void if the product is ever used in a commercial or business environment. The maximum liability of the warrantor under this warranty is limited to the purchase price actually paid by the customer for the product covered by the warranty, as confirmed by proof of purchase, regardless of the amount of any other direct or indirect damage suffered by the customer.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adaptation to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED IN DURATION TO THE DURATION OF THIS APPLICABLE WRITTEN WARRANTY. Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

IN NO EVENT SHALL THE WARRANTOR BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS, IMPLIED OR ANY OTHER THEORY OF LIABILITY, WHATSOEVER. Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you.

The warrantor does not authorize anyone, including, but not limited to, retailers, the subsequent consumer purchaser of the product from a retailer or remote purchaser, to obligate the warrantor in any way beyond the terms set forth herein.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resealed products, including but not limited to sale of such products on Internet auction sites and/or products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate in connection with any products or parts thereof which are repaired, replaced, altered, or modified, without the prior explicitly written consent of the warrantor.

This warranty gives you specific legal rights, and you may also have other rights which may vary from state to state.

Made in China.

Distributed by: CVS Pharmacy, Inc.

One CVS Drive, Woonsocket, RI 02895

© 2011 CVS/pharmacy

www.cvs.com 1-800-shop-CVS

For product information, questions or customer service, call toll free 1-800-536-0366.

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Distribuido por: CVS Pharmacy, Inc.